#### Class core values

- 1. Be **respect**ful to yourself and others
- 2. Be **confident** and believe in yourself
- 3. Always do your **best**
- 4. Be **cooperative**
- 5. Be **creative**
- 6. Have **fun**
- 7. Be **patient** with yourself while you learn
- 8. Don't be shy to **ask "stupid" questions**






#### Learning Objectives

- 1. Describe the concepts of selection and screening
- 2. Identify the main methods for evolving binders
- Critically evaluate the use of screening/selection methods for a given application
- 4. Identify methods for linking phenotype and genotype




#### The overall process of directed evolution





A key part of directed evolution is choosing the

best variants






**Screening** is the process of going through all variants and picking up the best one




#### Fluorescence can be used for screening






#### Fluorescence can be used for screening





#### Many enzymatic assays are screening-based






When you perform **selection**, only the fit variants survive



### Antibiotic plates are the most common method for selection






The outcome of evolution studies heavily depend on your selection/screening method




# The outcome of evolution studies heavily depend on your selection/screening method






### The outcome of evolution studies heavily depend on your selection/screening method






### In-class activity: You get what you select for ...





The selected variants need to be amplified to

create the next parents





### Selection happens at the protein level



# **Selection** happens at the protein level, but **amplification** is at the DNA level




## Selection happens at the protein level, but amplification is at the DNA level

We need a way to link the phenotype (protein) to genotype (DNA)




### In low-throughput assays ...





### In low-throughput assays, you can take note of the phenotype-genotype link



| Row ID | Activity | Sequence        |
|--------|----------|-----------------|
| A1     | 90%      | LKNMGFTHILKDFSA |
| A2     | 25%      | LKQMGFSHILKDWSA |
| A3     | 40%      | IRNMGYTHIVKDFSA |
|        |          |                 |
| H12    | 35%      | LRNCGWTHIIKDFTV |



**Containment**: Keeping the activity inside the cell where the DNA is!




### Antibiotic resistance is an example of contained activities






# Fluorescent protein signals are contained within the cytosol






## Oil-water emulsions can be used to contain the phenotype and genotype within the same droplet






**Display**: Covalently linking the protein to the DNA source












































# Yeast surface display is the method of choice for most protein engineering approaches







Yeast surface display can be linked with FACS to obtain best binders


































Proteins can be displayed on the surface of bacteria




## Proteins can be displayed on the surface of bacteria






Ribosome display can also be used for displaying

proteins



# mRNA display is attracting attentions as a method for selecting peptide binders





### Each method has its limitations/applications

|   | Phage display   | Eukaryotic<br>display         | Prokaryotic display                              | Ribosome<br>display | mRNA/cDNA<br>display |
|---|---|-------------------------------|--|---------------------|----------------------|
| Host organism   | Filamentous<br>phages,<br>M13, T4, T7,<br>lambda,<br>phagemid | S. cerevisiae,<br>P. pastoris | E. coli, B. subtilis,<br>L. bacillus, S. camosus | In vitro            | In vitro             |
| Library size Highest affinity $K_d$ (M) $^a$ Typical enrichment factor per round Nucleic acid selected Transformation required Library form |   |                               |  |                     |                      |
| Proteins to<br>be displayed   |   |                               |  |                     |                      |
| Covalent link<br>Surface anchorage  |   |                               |  |                     |                      |
| Post translational machinery  |   |                               |  |                     |                      |




### Each method has its limitations/applications

|  | Phage display  | Eukaryotic<br>display   | Prokaryotic display   | Ribosome<br>display  | mRNA/cDNA<br>display  |
|--|--|---|---|--|---|
| Host organism  | Filamentous<br>phages,<br>M13, T4, T7,<br>lambaa,              | S. cerevisiae,<br>P. pastoris   | E. coli, B. subtilis,<br>L. bacillus, S. camosus                            | In vitro   | In vitro  |
| T.:  | phagemid<br>10 <sup>9a</sup>                                   | 10 <sup>7</sup>   | 108-10  | 10 <sup>13-14</sup>  | 10 <sup>13-14</sup>   |
| Library size<br>Highest affinity<br>$K_d$ (M) <sup>a</sup> | 10<br>10 <sup>-11</sup>  | $10^{-14}$  | $10^{-13}$  | $10^{-12}$   | $10^{-10}$  |
| Typical enrichment factor per round                        | $10^{2-4}$   | 10 <sup>2-3</sup>   | $10^{2-3}$  | 10 <sup>1-3</sup>  | 10 <sup>1-3</sup>   |
| Nucleic acid<br>selected                                   | DNA  | DNA   | DNA   | mRNA   | mRNA/cDNA   |
| Transformation required                                    | Yes  | Yes   | Yes   | No   | No  |
| Library form   | Plasmid  | Plasmid   | Plasmid   | PCR fragment or mRNA   | mRNA/cDNA,<br>plasmid   |
| Proteins to<br>be displayed                                | Soluble, nontoxic,<br>compatible with<br>crossing<br>membranes | Soluble and<br>membrane,<br>nontoxic,<br>compatible<br>with crossing<br>membranes | Soluble and membrane,<br>nontoxic, compatible<br>with crossing<br>membranes | Most proteins including cytotoxic, chemically modified and membrane proteins | Soluble,<br>including<br>cytotoxic,<br>chemically<br>modified |
| Covalent link  | No   | No  | No  | No   | Yes (synthetic)   |
| Surface anchorage  | Capsid proteins  | Agglutination proteins, flocculation proteins                                     | Lpp-OmpA,<br>autotransporter<br>proteins, ice<br>nucleation<br>proteins     | Ribosome   | In vitro  |
| Post translational machinery                               | Simple   | Sophisticated   | Moderate  | Moderate   | Simple  |




# Continuous evolution methods combine diversification, selection and amplification





### Phage-Assisted Continuous Evolution (PACE)





#### For the next lecture:

- Pre-class assessment for the next lecture
   Needs to be done before the start of class, will be available after this class
- Post-class assignment Write up questions for our panelists
- 3. Second journal: Will be discussed next week



# Next lecture: The challenging case of enzymes

